

Lovibond® Water Testing Tintometer® Group



Water Safety Kit Microbiology Duo

Mobile microbiological test equipment with 2 incubators



Disaster situation or potable water monitoring in demanding situations: You need certainty and reliability for test results in everywhere. The Water Safety Kit Microbiology Duo contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria - with both incubators simultaneously! Perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. With the optional backpack you can take the kit everywhere with you.

- Testing coliforms and E. coli bacteria according to WHO standards
- Testing different forms of bacteria simultaneously with 2 incubators
- DI 20 with heating and cooling as standard

Ürün No: 56K681255

Sanayi

Diğer Sanayiler | Sivil Toplum Kuruluşları

Uygulama

İçme Suyu Arıtma

Water Safety Kit Microbiology Duo

In one hand or on the back: This test kit accompanies you in every situation with the complete microbiological test equipment. Test different sorts of bacteria directly on-site - with both incubators even simultaneously! Be sure of the quality of potable water everywhere, even when there is no laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

Teknik Veri

Boyutlar	487 x 229 x 386 mm
Ambalajlı Ağırlık	18.75 kg

Teslimat Kapsam Metni

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500ml
- Tweezers
- Pump, 35ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- 2 DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin)
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

Aksesuarlar

Başlık	Ürün No
Güç kablosu seti	136300
Incubator Battery, Lead/Acid	192345
Şişeler	363590
Universal container with spoon and lid, 30 ml	366239
Plastik kap, 100 ml derecelendirmeli, kapaklı	385130
Petri kapları, 10 adet.	400855
Su Güvenlik Kitleri için sırt çantası	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Ölçü kabı	56A018501
Diyaframlı filtreleme	56A023801
Cımbız	56A604201
Diyaframlı filtreleme	56B000710
D004 TTC/MAC daldırma kızakları (Macconkey)	56B010410
D005 TTC/Kromojenik daldırma kızakları	56B010510
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Almanya

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Birleşik Krallık

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Çin

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malezya

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brezilya

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Amerika Birleşik Devletleri

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Hindistan

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
İspanya