

Lovibond® Water Testing Tintometer® Group



Sewage / Effluent Test Kit

Test all important parameters for sewage/effluent analysis



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Ürün No: 56K701220

Sanayi

Diğer Sanayiler

Uygulama

Dezenfeksiyon Kontrol | Soğutma Suları

Teslimat Kapsam Metni

- Plastic case
- Checkit Comparator
- 4 - 10 pH Universal Checkit Disc
- Universal pH tablets
- Permanganat tablets
- Acidifying Sewage Effluent tablets
- 2 plastic vials 10 mL with lid
- Shaker tube with lid 100 mL
- Syringe plastic, 20 mL
- Titration jar with cap, plastic, 60 mL
- Turbidity tube
- Instruction manual

Aksesuarlar

Başlık	Ürün No
CHECKIT® diski pH (üniversal)	146130
Plastik kap, 100 ml derecelendirmeli, kapaklı	385130A
Kapaklı titrasyon kabı, plastik, 60 mL	56A006701
Bulanıklık testi 5-500 JTU	56A010801
Turbidity tube (one piece) - 30-400 NTU	56A607101

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Almanya

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Birleşik Krallık

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Çin

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malezya

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brezilya

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Amerika Birleşik Devletleri

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Hindistan

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
İspanya