

# Lovibond® Water Testing Tintometer® Group



## Water Treatment Engineers Test Kit (advanced)

Advanced test equipment for water treatment professionals



- All instruments & reagents for daily routine
- Including advanced photometer MD640 for additional parameters
- Early prevention & detection of risks
- Ensure system performance & operating costs

Ürün No: 56K701500

Water treatment engineers are engaged daily in the control of water quality in closed circuits, evaporative cooling systems and boiler systems. They need reliable equipment and accurate results to rely on. Both have an important impact on decisions on treatment programs and optimal system performance

This Advanced Engineers Kit contains all you need to protect the life of the plant and operating costs. It allows early detection of corrosion, scale and biofouling and the investigation of increased biological activity.

Includes the following parameters: Conductivity, pH, glycol, aluminum, bromine, chlorine, copper, iron, molybdate, phosphate and alkalinity (M, P, OH).

### Sanayi

Belediyeler | Deniz Sanayi | Diğer Sanayiler | Güç ve Enerji | Halka-Açık Havuz | Havuz Hizmetleri | İlaç Sanayi | Kimya Sanayi | Sivil-Toplum Kuruluşları | Yağ ve Gaz | Yerel Havuz | Yiyecek ve İçecek-Sanayi

### Uygulama

Atık Su Arıtma | Dezenfeksiyon Kontrol | Galvanizasyon | Ham Su-Arıtma | Havuz Suyu Kontrol | İçme Suyu Arıtma | Kazan Suları | Ot-hers | Soğutma Suları

### Water Treatment Engineers Test Kit (advanced)

This kit contains pocket meters, checkit comparator, drop tests and MD640 photometer for advanced testing an extended range of control parameters.

Take the complete equipment for direct on-site analysis in practical case in one hand.

Advances analysis of all types of industrial water with visual, electrochemical and photometric methods.

## Ölçüm Aralığı

Test Name	Ölçüm Aralığı	Kimyasal Metod
Alkalinity (P, M, OH)	50 - 2400 mg/L CaCO <sub>3</sub>	
Alüminyum T	0.01 - 0.3 mg/L Al	Eriokrom Cyanine R
Bakır T	0.05 - 5 mg/L Cu <sup>®</sup>	Biquinoline
Brom T	0.05 - 13 mg/L Br <sub>2</sub>	DPD
Chloride	20 - 12000 mg/L Cl <sup>-</sup>	
Conductivity	0 - 200 mS/cm	
Demir HR L	0.1 - 10 mg/L Fe	Tiyoglikolat
Demir LR L (A)	0.03 - 2 mg/L Fe	Ferrozin/Tiyoglikolat
Fosfat HR T	0.33 - 26 mg/L P	Vanadomolybdate
Glikol	% PEG/MEG	
Hardness Calcium	5 - 600 mg/L CaCO <sub>3</sub>	
Hydrogen Peroxide	15 - 500 mg/L H <sub>2</sub> O <sub>2</sub>	
Klor T	0.01 - 6.0 mg/L Cl <sub>2</sub> <sup>®</sup>	DPD
Molibdat T	1 - 50 mg/L MoO <sub>4</sub>	Tiyoglikolat
Nitrit	10 - 2000 mg/L NaNO <sub>2</sub>	
pH	-2 - 16	
Phosponate	0 - 20 mg/L HEDP	
PTSA	10 - 400 ppb	Floresan
Sertlik (Evet/Hayır)	8 - 20 mg/L CaCO <sub>3</sub>	
Sulphite	25 - 150 mg/L Na <sub>2</sub> SO <sub>3</sub>	
Tannin	50 - 300 mg/L Tannin	
Toplam sertlik derecesi	5 - 600 mg/L CaCO <sub>3</sub>	

## Teslimat Kapsam Metni

- Plastic case
- 6 Cuvettes 24 mm round with lids & sealing rings
- SD 335 Multi
- MD640 Photometer
- Brush, 11 cm length
- DPD No.1 tablets BT
- DPD No.3 tablets
- Glycine tablets
- DPD Nitrite tablets bottle
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No.1 Test tablets
- Copper No.2 Test tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- Phosphate H.R.P1 tablets
- Phosphate H.R.P2 tablets
- Copper/Zinc L.R.Test Tablets
- Hardness Yes/No tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- 1413uS Standard Conductivity Solution
- Iron Reagent FE6
- Iron Reagent FE5
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Alkalinity OH Reagent
- Hardness Total Buffer TH2
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- Tannin Indicator TN1
- Tannin Titrant TN2
- Acetic acid 20%
- Nitrite Titrant N2
- Phosponate Titrant P5
- Hydrogen Peroxide Buffer HP1
- Phosponate Indicator P4L
- Nitrite Indicator N1
- Hydrogen Peroxide HR Titrant HP2
- Iron Reagent FE8
- Phosponate pH Adjuster P3
- Chloride Indicator CC1
- Calcium Hardness Indicator CH1P
- Calcium Hardness Buffer CH2
- Phosponate Neutraliser P1/2
- Sulphite Indicator S1
- Sulphite Indicator S2

## Aksesuarlar

Başlık	Ürün No
SD 320 / SD 325 Con için LC 12, 200 mS/cm'ye kadar	19805040
ValidCheck Klor 1,5 mg/l	48105510
Enjektör, plastik, 20 mL	56A006501
Kapaklı titrasyon kabı, plastik, 60 mL	56A006701
Plastik enjektör, 5 mL	56A008501
Plastik enjektör, 1 mL	56A013501
Test Tüpü 5/10 mL + CAP	56A600401
Glikol refraktometresi	56M000401
SD 335 Multi (ana cihaz)	724830

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Almanya

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Birleşik Krallık

### Tintometer China

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Çin

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malezya

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brezilya

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Amerika Birleşik Devletleri

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Hindistan

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
İspanya