

Lovibond® Water Testing

Tintometer® Group



5in1 Mini

The focus on various parameters simultaneously



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 151700

Inexpensive

The Pool Test Kit is the cost-effective entry-level variant in water analysis.

Variable

The different versions of the Pool Test Kit allow you to test different parameters. Possible combinations are, for example, chlorine + pH, bromine + pH, active oxygen + pH or biguanide (PHMB) + pH.

Ready for immediate use

The Pool Test Kit includes the tablet reagents for 20 tests each and instruction manual, all packaged in a handy and transparent mini suitcase. So you not only keep track of everything, but also have the test kit at your fingertips at any time.

Industry

Pool Domestic | Pool Services

Application

Disaster Management | Pool Water Control

5in1 Mini

Perfect for the entry into the water analysis in the private environment. Measure several parameters at the same time.

Measuring Range

Test Name	Measuring Range
Alkalinity-m	50 - 300 mg/L CaCO ₃
Bromine	0.2 - 6.8 mg/L Br ₂
Chlorine	0.1 - 3.0 mg/L Cl ₂
Cyanuric acid	20 - 200 mg/L CyA
pH	6.8 - 8.2 pH

Delivery Scope

- Pool Test Kit in a bubble pack or mini case
- Instruction manual

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain