

Lovibond® Water Testing

Tintometer® Group



Trousse d'analyse à l'isothiazolinone, de 0 à 7.5 mg/L

For the control of Isothiazoline concentration



- For monitoring of cooling water & other industrial applications
- Works with colour card
- For up to 100 tests
- Carry out the test in a few steps directly on site

Code: 56K701200

Biocide with versatile effect

Isothiazolinone is a biocide for the control of bacteria, fungi and algae, especially in industrial water and cooling water. It works by separating the proteins even at low doses. This inhibits microbial growth and kills the unwanted "inhabitants". Isothiazolinone has good degradability and no side effects. It is safe to use and has good stabilisation properties. Its use is inexpensive.

Easy to use

Professionals and non-professionals can quickly obtain reliable results with this Test Kit. The test works with the colourimetric method: If isothiazolinone is present in the water, the sample turns blue after the addition of certain reagents. The concentration can be read off from the colour chart. The test can be carried out directly on site and provides certainty within minutes - without a laboratory. The easy-to-understand instructions with pictograms ensure that the test can be carried out by any user, regardless of experience.

Many advantages

Many advantages With our Test Kit you can quickly determine how much active isothiazolinone is effective in the water. This helps you avoid incorrect dosing and possible material damage or system blockages.

Industrie

Autres industries

Applications

Contrôle de la désinfection | Eau de refroidissement

Trousse d'analyse à l'isothiazolinone, de 0 à 7.5 mg/L

Isothiazolinone has a broad spectrum of efficacy against microbiological growth in industrial process water, especially in cooling water. It is already effective at very low concentrations. To ensure effectiveness, regular checks of the isothiazolinone concentration are recommended - to avoid costs and possible damage to the material.

The Lovibond® Non-Oxidising Biocide Test Kits also provide a simple solution to ensure bleed concentrations are in line with environmental consents.

Gamme de mesure

Test Name	Gamme de mesure
Isothiazoline	0 - 7.5 mg/L

Contenu de livraison

- Plastic case
- Isothiazolinone DK1
- Isothiazolinone DK2
- Isothiazolinone DK3
- Isothiazolinone DK4
- Isothiazolinone DK5
- Stirring rod and powder spoon
- Plastic tube with lid, 10 mL
- Syringe plastic, 20 mL
- 2 Titration jars with cap, plastic, 60 mL
- Instruction manual

Accessoires

Titre	Code
Seringue, plastique, 20 mL	56A006501
Agitateur et cuiller à poudre	56A006601
Récipient de titrage avec couvercle, plastique, 60 mL	56A006701
Tube en plastique avec couvercle 10 ml	56A006901
Tube de test 5/10 mL + capuchon	56A600401

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Allemagne

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Royaume-Uni

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Chine

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaisie

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brésil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
États-Unis

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Inde

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Espagne