

# Lovibond® Water Testing Tintometer® Group



## Water Safety Kit Combined

Special test kit for potable water not only for NGOs



- Emergency response kits for disaster situations and water monitoring outback
- Equipment for chemical and microbiological analysis in two cases
- Also available as a duo version

货号: 56K681251

You need safe, fast and reliable test results in difficult situations like disasters or potable water monitoring in rural environments without laboratories? You immediate certainty about the quality of potable water to save health and lives? Take all you need with you and combine microbiological methods and simple chemical methods of analysis.

Take the smaller Water Safety Kit Chemical in one hand and the bigger Water Safety Kit Microbiology on your back in the optional backpack. You are able to test different sorts of bacteria on-site due to the heating and cooling function of our DI 20 petri dish incubators.

### 行业

NGO | 公共泳池 | 其他行业 | 化工行业 | 医药行业 | 家庭泳池 | 市政 | 油气 | 泳池服务 | 海洋工业 | 能源 | 食品和饮料行业

### 应用

冷却水 | 原水处理 | 污水处理 | 泳池水质控制 | 消毒控制 | 锅炉水 | 饮用水处理

### Water Safety Kit Combined

Combine both microbiological and physio chemical analysis of some of the most critical parameters to determine the suitability of water for drinking. Take the whole equipment with you everywhere you need immediate and reliable test result.

With their rugged, waterproof cases the kits are easy to transport. Even non technical personal is able to use it due to pictorial instructions.

### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

## 测量范围

Test Name	测量范围	化学方法
Chlorine T	0.01 - 6.0 mg/L Cl <sub>2</sub> a)	DPD
Conductivity	0 - 20 mS/cm	
pH	0 - 14	
T 硝酸盐	0.08 - 1 mg/L N	锌还原 / NED
T 铵	0.02 - 1 mg/L N	靛酚蓝
好氧菌 (TVC) / 大肠杆菌、大肠菌群		
好氧菌 (TVC) / 大肠菌群		
浊度	0.01 - 1100 NTU	

## 标准配置文本

- 1 case Water Safety Kit Microbiology
- 1 case Water Safety Kit Chemical

标题	货号
Media Pads and Dispenser (Pack of 200)	56A016330
量杯	56A018501
膜式过滤	56A023801
镊子	56A604201
膜式过滤	56B000710
D004 TTC/MAC Dipslide (麦康凯)	56B010410
D005 TTC/呈色法 Dipslide	56B010510
KS3-pH 10 缓冲溶液 65 ml	56L000365
KS6-pH4 缓冲溶液 65 ml	56L000665
KS9-pH 7 缓冲溶液 65 ml	56L000965
KS16-1413US 标准电导率溶液	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

## 辅料

标题	货号
电源线组	136300
Incubator Battery, Lead/Acid	192345
浊度标准 T-CAL 套件 (<0.1, 20, 200, 800 NTU)	194150
瓶子	363590
Universal container with spoon and lid, 30 ml	366239
塑料容器, 100 ml 刻度, 带盖	385130
塑料容器, 30 ml 刻度, 带盖	385132
培养皿, 10个。	400855
ValidCheck 氯 1.5 mg/l	48105510
搅拌棒和粉勺	56A006601
水上安全套装背包	56A01400

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
德国

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
英国

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
中国

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebu Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
马来西亚

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
巴西

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
美国

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
印度

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
西班牙