

# Lovibond® Water Testing Tintometer® Group



## Water Safety Kit Microbiology Duo

Mobile microbiological test equipment with 2 incubators



Disaster situation or potable water monitoring in demanding situations: You need certainty and reliability for test results in everywhere. The Water Safety Kit Microbiology Duo contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria - with both incubators simultaneously! Perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. With the optional backpack you can take the kit everywhere with you.

- Testing coliforms and E. coli bacteria according to WHO standards
- Testing different forms of bacteria simultaneously with 2 incubators
- DI 20 with heating and cooling as standard

货号: 56K681255

### 行业

NGO | 其他行业

### 应用

饮用水处理

### Water Safety Kit Microbiology Duo

In one hand or on the back: This test kit accompanies you in every situation with the complete microbiological test equipment. Test different sorts of bacteria directly on-site - with both incubators even simultaneously! Be sure of the quality of potable water everywhere, even when there is no laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

### Technical Data

|        |                    |
|--------|--------------------|
| 尺寸     | 487 x 229 x 386 mm |
| 带包装的重量 | 18.75 kg           |

## 标准配置文本

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500ml
- Tweezers
- Pump, 35ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- 2 DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

## 辅料

| 标题  | 货号        |
|---|-----------|
| 电源线组  | 136300    |
| Incubator Battery, Lead/Acid                  | 192345    |
| 瓶子  | 363590    |
| Universal container with spoon and lid, 30 ml | 366239    |
| 塑料容器, 100 ml 刻度, 带盖                           | 385130    |
| 培养皿, 10个。                                     | 400855    |
| 水上安全套装背包                                      | 56A01400  |
| Media Pads and Dispenser (Pack of 200)        | 56A016330 |
| 量杯  | 56A018501 |
| 膜式过滤  | 56A023801 |
| 镊子  | 56A604201 |
| 膜式过滤  | 56B000710 |
| D004 TTC/MAC Dipslide ( 麦康凯 )                 | 56B010410 |
| D005 TTC/呈色法 Dipslide                         | 56B010510 |
| Membrane Lauryl Sulphate Broth, 38.1 gr.      | 56P069740 |
| Pump (35 ml stroke)                           | 56X610301 |
| Pipette, disposable, 1.5 ml                   | BI-K10060 |
| Marker Pen                                    | BI-K14087 |

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
德国

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
英国

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
中国

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
马来西亚

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
巴西

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
美国

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
印度

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
西班牙