

Lovibond® Water Testing Tintometer® Group



Sewage / Effluent Test Kit

Test all important parameters for sewage/effluent analysis



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

货号: 56K701220

行业

其他行业

应用

冷却水 | 消毒控制

辅料

标题	货号
CHECKIT® 色盘 pH (常规)	146130
塑料容器, 100 ml 刻度, 带盖	385130A
带盖滴定容器, 塑料, 60毫升	56A006701
浊度测试 5-500 JTU	56A010801
Turbidity tube (one piece) - 30-400 NTU	56A607101

标准配置文本

- Plastic case
- Checkit Comparator
- 4 - 10 pH Universal Checkit Disc
- Universal pH tablets
- Permanganat tablets
- Acidifying Sewage Effluent tablets
- 2 plastic vials 10 mL with lid
- Shaker tube with lid 100 mL
- Syringe plastic, 20 mL
- Titration jar with cap, plastic, 60 mL
- Turbidity tube
- Instruction manual