

Lovibond® Water Testing Tintometer® Group



Water Treatment Engineers Test Kit (visual)

Monitoring water treatment for professionals with visual methods



The Lovibond® Visual Engineers Test Kit is designed for the water treatment professional to be used on all types of industrial water systems such as cooling towers, closed water systems and steam raising boilers.

The Kit uses visual methods on a comparator. It includes drop tests and pocket meters for electrochemistry parameters.

Each test has been chosen specifically for these applications to ensure accurate repeatable results.

Includes the following parameters: Conductivity, pH, glycol, aluminium, bromine, chlorine, copper, iron, molybdate, phosphate.

- Mobile laboratory for water treatment monitoring
- Works with visual tests & further methods
- Ensure system performance & operating costs
- Most efficient monitoring of industrial water samples

货号: 56K701300

行业

其他行业 | 化工行业 | 医药行业 | 市政 | 油气 | 海洋工业 | 能源 | 食品和饮料行业

应用

Others | 冷却水 | 污水处理 | 消毒控制 | 锅炉水

Water Treatment Engineers Test Kit (visual)

You need safe, fast and reliable test results?

Take everything you need for chemical analysis of industrial water quality in one hand

This comprehensive Test Kit combines everything water professionals need for on-site analysis with visual test methods and more.

测量范围

Test Name	测量范围	化学方法
Chloride	20 - 12000 mg/L Cl ⁻	
Conductivity	0 - 20 mS/cm	
Glycol	% PEG/MEG	
Hardness Calcium	5 - 600 mg/L CaCO ₃	
Hydrogen Peroxide	15 - 500 mg/L H ₂ O ₂	
Iron LR T	0.05 - 1 mg/L Fe	PPST
pH	0 - 14	
T 铝	0 - 0.3 mg/L Al	依来铬靛蓝 R
亚硝酸盐	10 - 2000 mg/L NaNO ₂	
亚硫酸盐	25 - 150 mg/L Na ₂ SO ₃	

Test Name	测量范围	化学方法
总硬度	5 - 600 mg/L CaCO ₃	
氯 MR T	0 - 2 mg/L Cl ₂	DPD
溴 T	0 - 5 mg/L Br ₂	DPD
硬度 (是/否)	8 - 20 mg/L CaCO ₃	
碱度 (P, M, OH)	50 - 2400 mg/L CaCO ₃	
磷酸盐 HR T	0 - 80 mg/L PO ₄	Vanadomolybdate
膦酸盐	0 - 20 mg/L HEDP	
钼酸盐 HR T	50 - 500 mg/L MoO ₄	巯基乙酸
铜 HR T	0 - 5 mg/L Cu	双喹啉
鞣	50 - 300 mg/L Tannin	

标准配置文本

- Plastic case
- SD40 Multi Pocket Tester
- Checkit Comparator
- 0 - 5 mg/L Bromine Checkit® Disc
- 0 - 4 mg/L Chlorine DPD Checkit® Disc
- 0 - 1.0 mg/L Iron LR Checkit® Disc
- 50 - 500 mg/L Molybdate HR Checkit® Disc
- 0 - 0.3 mg/L Aluminium Checkit® Disc
- 0 - 5.0 mg/L Copper HR Free & Total Checkit® Disc
- 0 - 80 mg/L Phosphate HR Checkit® Disc
- 1 - 10 mg/L Iron HR Checkit® Disc
- 2 Plastic vials 10 mL with lid
- DPD No.1 Rapid tablets BT
- DPD No.3 Rapid tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No.1 Test tablets
- Copper No.2 Test tablets
- Phosphate HR Comp. tablets
- Hardness Yes/No tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- 1413uS Standard Conductivity solution
- Iron Reagent FE5
- Iron Reagent FE6
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Acidity / Alkalinity P Indicator PA1
- Alkalinity OH Reagent
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- Tannin Indicator TN1
- Hydrogen Peroxide HR Titrant HP2
- Acetic acid 20%
- Nitrite Titrant N2
- QAC Buffer QA2
- Phosphonate Indicator P4L
- Nitrite Indicator N1
- Tannin Titrant TN2
- Iron Reagent FE8
- Phosphonate pH Adjuster P3
- Chloride Indicator CC1
- Calcium Hardness Indicator CH1P
- Calcium Hardness Buffer CH2

辅料

标题	货号
CHECKIT® 比色仪	145000
塑料注射器 20 mL	56A006501
带盖滴定容器, 塑料, 60毫升	56A006701
塑料注射器, 5 mL	56A008501
塑料注射器 1 ml	56A013501
试管 5/10 mL + CAP	56A600401
乙二醇折光仪	56M000401

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
德国

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
英国

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
中国

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
马来西亚

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
巴西

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
美国

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
印度

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
西班牙