


Lovibond® Water Testing

Tintometer® Group



Closed System Engineers Test Kit

The expert tool for reliable result



These Test Kits for cooling water monitoring are specially designed for engineers and mobile use. It provides the equipment and reagents required to monitor the correct amount and dosing levels of various important parameters. You can carry everything you need with you and test directly onsite. The Kit is provided in a robust carry case and includes all important equipment from measurement instruments to reagents.

- Complete analysis kit including molybdate, iron and conductivity
- Prevents corrosion and scale
- Accurate inhibitor dosage rates
- Reproducible, accurate results

货号: 56K701600

行业

其他行业 | 化工行业 | 医药行业 | 市政 | 油气 | 海洋工业 | 能源 | 食品和饮料行业

应用

冷却水 | 污水处理 | 锅炉水

Closed System Engineers Test Kit

Due to its heat absorption capacity, water is the preferred cooling medium in many areas of industry. However, water always contains many dissolved and undissolved substances that can cause trouble even when travelling in closed cooling circuits. Extensive water treatment or chemical post-treatment by experts is always necessary. Problems are mainly caused by deposits, corrosion and microbiological activities. For the experts, a regular check of certain parameters is therefore just as essential as reliable equipment.

测量范围

Test Name	测量范围	化学方法
Alkalinity M	50 - 2400 mg/L CaCO ₃	
Chloride	20 - 12000 mg/L Cl ⁻	
Glycol	% PEG/MEG	
Iron LR T	0.05 - 1 mg/L Fe	PPST
T 铝	0 - 0.3 mg/L Al	依来铬氰蓝 R
亚硝酸盐	10 - 2000 mg/L NaNO ₂	
总硬度	5 - 600 mg/L CaCO ₃	
钼酸盐 HR T	50 - 500 mg/L MoO ₄	巯基乙酸
铜 HR T	0 - 5 mg/L Cu	双喹啉

标准配置文本

- Plastic case
- SD40 Multi Pocket Tester
- Checkit Comparator
- 50 - 500 mg/L Molybdate HR+ Checkit Disc
- 0 - 1.0 mg/L Iron LR Checkit Disc
- 0 - 0.3 mg/L Aluminium Checkit Disc
- 0 - 5.0 mg/L Copper HR Free and Total Checkit Disc
- 2 Plastic vials 10 mL with lid
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No. 1 Test tablets
- Copper No. 2 Test tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- Iron Reagent FE5
- Hardness Total Buffer TH2
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- 1413µS Standard Conductivity Solution
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Nitrite Indicator N1
- Nitrite Titrant N2
- Chloride Indicator CC1
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- 2 Titration jars with caps, plastic, 60 mL
- Syringe plastic, 20 mL
- Wash Bottle, 250 mL
- Syringe Filter Holder 25 mm
- Filter Circle 0.45 µm, 25 mm
- Glycol Refractometer, RHA-218ATC
- Instruction manual

辅料

标题	货号
CHECKIT® 比色仪	145000
塑料注射器 20 mL	56A006501
带盖滴定容器, 塑料, 60毫升	56A006701
塑料注射器, 5 mL	56A008501
洗瓶 250 ml	56A009001
塑料注射器 1 ml	56A013501
圆滤光片 0.45 um , 25 mm	56A020050
试管 5/10 mL + CAP	56A600401
乙二醇折光仪	56M000401

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
德国

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
英国

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
中国

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
马来西亚

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
巴西

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
美国

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
印度

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
西班牙