

Lovibond® Water Testing

Tintometer® Group



废水测量站 MD 610

Waste water control with mobile data sets



- Mobile photometer MD 610
- Memory for 500 measurement results
- Data transmission via **Bluetooth®** and **App AquaLX®**
- Thermoreactor RD 125 with 24 tube holes

货号: 214110

行业

NGO | 其他行业 | 化工行业 | 医药行业 | 市政 | 油气 | 海洋工业 | 能源 | 食品和饮料行业

应用

原水处理 | 污水处理 | 电镀 | 饮用水处理

废水测量站 MD 610

You have to react immediately to the daily control routine - no matter where and when? With the mobile photometer MD 610, you can call up, evaluate and share the measurement results of the wastewater directly on site on your smartphone or tablet with your team. The interference filters ensure increased reproducibility of the results. With 120 pre-programmed methods and memory space for around 500 measured values, you are flexible in your routine work. With the RD 125 thermoreactor with 24 holes for 16 mm tubes, precisely fitting reagents for COD measurement and the appropriate accessories, you get the complete waste water set-up from a single source and save costs.

测量范围

Test Name	测量范围	化学方法
HR TT 铰	1.0 - 50 mg/L N	水杨酸
LMR TT 化学需氧量	15 - 300 mg/L COD ^{b)}	Dichromate / H ₂ SO ₄
LR TT 化学需氧量	3 - 150 mg/L COD ^{b)}	Dichromate / H ₂ SO ₄
MR TT 化学需氧量	20 - 1500 mg/L COD ^{b)}	Dichromate / H ₂ SO ₄
PP 亚硝酸盐	0.01 - 0.3 mg/L N	重氮化
TN HR TT	5 - 150 mg/L N ^{b)}	过硫酸盐消解
TT 硝酸盐	1 - 30 mg/L N	变色酸
TT 磷酸盐水合物	0.02 - 1.6 mg/L P ^{b)}	磷钼蓝

Technical Data

便携性 便携式

标准配置文本

- 光度计配有标配配件
- 热反应器 RD 125
- COD 比色杯测试 0 - 150 mg/l 和 0 - 1500 mg/l
- 铵 VARIO HR 比色杯测试 1 - 50 mg/l N
- 硝酸盐 VARIO 比色杯测试 1 - 30 mg/l N
- 亚硝酸盐 LR VARIO 粉包 0.01 - 0.3 mg/l N
- 氮 VARIO 总 HR 比色杯测试 5 - 150 mg/l N
- 磷酸盐 VARIO 总 HR 比色杯测试 0.02 - 1 mg/l P/0.06 - 3.5 mg/l PO₄
- 比色杯架
- 膜式过滤组
- 操作说明
- 保修声明

辅料

标题	货号
带盖的圆形比色杯 Ø24mm, 高 48 mm, 10 ml, 12 件套	197620
圆形比色杯, 盖子直径 Ø 24 mm, 高 48 mm, 10 ml, 5 件套	197629
圆形比色杯, 盖子直径 Ø 16 mm, 高 90 mm, 10 ml, 10 件套	197665
圆形比色杯 16 mm 适配器	19802190
圆形比色杯 13 mm 适配器	19802192
热反应器 RD 125	2418940
用于样本处理的薄膜过滤套件, 25 个薄膜过滤器 0.45 微米, 2 个注射器 20 ml	366150
试样瓶支架, 用于 6 个 Ø 24 mm 圆形试样瓶	418951
用于 10 个 Ø 16 mm 圆形比色杯的比色杯支架	418957
移液器吸头, 1-5 ml (白色) 100 件	419066
移液器吸头, 0.1-1 ml (蓝色), 1000 件	419073
自动移液器, 1-5 ml	419076
自动移液器, 0.1-1 ml	419077
带把手的塑料漏斗	471007
W100/OG/10MM 矩形比色杯, 光学玻璃	601040
W100/OG/50MM 矩形比色杯, 光学玻璃	601070
W110/UV/10MM 矩形比色杯, 石英紫外线	661130
Factory calibration certificate ISO 9001 for MD600/MD610/MD640/MultiDirect	999752

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
德国

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
英国

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
中国

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
马来西亚

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
巴西

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
美国

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
印度

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
西班牙